

TOP CUT

TRU-CUT NEEDLE FOR AUTOMATIC DEVICES

The **TOP CUT** is a tru-cut biopsy needle for **PRO-MAG™** type automatic systems.

The **TOP CUT** needles features are:

- Marked needle with Menghini type tru-cut point
- Whistle pointed stylet
- Inner echogenic marker
- 16cm notch
- 11-22mm advancement
- It includes a spacer that makes easier the needle positioning inside the automatic system
- Colour code

TOP CUT - 22mm advancement

CODE	Diameter	Length	Box
TOP 1410	14G	100mm	10
TOP 1415	14G	150mm	10
TOP 1420	14G	200mm	10
TOP 1610	16G	100mm	10
TOP 1615	16G	150mm	10
TOP 1620	16G	200mm	10
TOP 1810	18G	100mm	10
TOP 1815	18G	150mm	10
TOP 1820	18G	200mm	10
TOP 1825	18G	250mm	10
TOP 1830	18G	300mm	10
TOP 2010	20G	100mm	10
TOP 2015	20G	150mm	10
TOP 2020	20G	200mm	10

TOP CUT Short 11mm advancement

CODE	Diameter	Length	Box
TOP 14075S	14G	75mm	10
TOP 14115S	14G	115mm	10
TOP 16075S	16G	75mm	10
TOP 16115S	16G	115mm	10
TOP 18075S	18G	75mm	10
TOP 18115S	18G	115mm	10

*PRO-MAG is an MD TECH's registered trade-mark

CE 0373



HISTOLOGICAL BIOPSY

H.S. Hospital Service S.p.A. is certified ISO 9001 and EN 46001



H.S. HOSPITAL SERVICE S.p.A.
Via delle Valli snc - 04011 Aprilia (LT)
Tel. 069201961 - Fax 0692727871 - Orders fax 069275519
www.hshospitalservice.com e-mail: hs@hshospitalservice.com