

UROCORE

TRU-CUT TYPE MANUAL DEVICE FOR HISTOLOGICAL BIOPSY

The **UROCORE** is a tru-cut manual needle for histological biopsy.

The depth marks allow a safe positioning of the needle.

The luer hub makes possible a cytological aspiration through a syringe after the histological biopsy.

The removable inner stylet features a 2cm notch that can be removed.

The inner echogenic marker prevents the spreading of cancerous cells.

CODE	Diameter	Length	Box
UC 14115	14G	115mm	10
UC 1415	14G	150mm	10
UC 1420	14G	200mm	10
UC 16115	16G	115mm	10
UC 1615	16G	150mm	10
UC 1620	16G	200mm	10
UC 18115	18G	115mm	10
UC 1815	18G	150mm	10
UC 1820	18G	200mm	10
UC 21115	21G	115mm	10
UC 2115	21G	150mm	10
UC 2120	21G	200mm	10



CE 0373

HISTOLOGICAL BIOPSIY

H.S. Hospital Service S.p.A. is certified ISO 9001 and EN 46001



H.S. HOSPITAL SERVICE S.p.A.

Via delle Valli snc - 04011 Aprilia (LT)
Tel. 069201961 - Fax 0692727871 - Orders fax 069275519
www.hshospitalservice.com e-mail: hs@hshospitalservice.com